

COVID-19 Antibody Testing

Anywhere, Anytime for \$30

➤ Are you tired of the slow COVID-19 testing process?

- Long wait times for an appointment
- Testing at offsite location
- Week long wait to obtain results, and quarantine during that time
- You, your employees, and/or students lose valuable time

Introducing the i-calQ Solution:



The i-calQ Solution brings diagnostic testing direct to you with results in **10 minutes**. Our technology uses smartphones to analyze and quantify antibody levels. This portable solution can be used by companies, schools, and at sporting events to triage patients on the spot to help reduce the spread of COVID-19.

➤ How it Works:

- The i-calQ Solution uses lateral flow immunoassay test cassettes with Emergency Use Authorization through the FDA
- A simple finger prick is all you need to draw two drops of blood to put into the cassette
- The cassette is then placed onto the i-calQ reader and our proprietary App analyzes and quantifies the antibody levels, just like centralized labs do, except our point of care device can be used **Anywhere, Anytime, Affordably with results in just 10 minutes.**

➤ Who should use this device:

- The i-calQ Solution is perfect for any organization, who wants to take control of COVID-19 testing.
- With this device you can test daily, weekly, or as often as you want to ensure you, your employees and students are receiving fast efficient test results to keep your organization safe.

With 7 issued patents in the USA, i-calQ is the only company that can provide this smartphone-based solution.

To learn more please visit our website at: <https://i-calQ.com>.

